

If digital image analysis is used, additional requirements in the Digital Image Analysis section also apply.

****REVISED** 08/24/2023**

CYG.47880 Report Elements

Phase I

For in situ hybridization (ISH) tests that provide independent predictive information, the patient report includes information on specimen processing, the probe, and the scoring method used.

NOTE: The following information must be included in the patient report:

1. The type of specimen fixation and processing (eg, formalin-fixed paraffin-embedded sections, air-dried imprints)
2. The probe and, if applicable, the detection system used (ie, LSAB, polymer, proprietary kit, vendor name, etc.; information on the type of equipment used is not necessary)
3. Criteria used to determine a positive vs. negative result and scoring system (eg, manual or automated)
4. Laboratory interpretation of predictive marker testing (ISH) is reported according to the manufacturer's instructions, or when available, following the structure, format, and criteria set forth in the current [CAP guidelines](#) relating to predictive marker testing (eg, ASCO/CAP HER2 testing in breast cancer and CAP/ASCP/ASCO HER2 in gastroesophageal carcinoma).
5. Limitations relating to suboptimal preanalytical factors that may impact results, such as prolonged cold ischemia time, unknown ischemia time, or over- or under-fixation.

Evidence of Compliance:

- ✓ Report template containing all required elements **AND**
- ✓ Copies of patient reports confirming inclusion of the required elements **AND**
- ✓ Established guidelines used by the laboratory

REFERENCES

- 1) Wolff AC, Somerfield MR, Dowsett M, et al. Human Epidermal Growth Factor Receptor 2 Testing in Breast Cancer: American Society of Clinical Oncology/College of American Pathologists Guideline Update. *Arch Pathol Lab Med*. Published online June 7, 2023. doi: 10.5858/arpa.2023-0905-SA.
- 2) Bartley AN, Washington MK, Ventura CB, et al. HER2 Testing and Clinical Decision Making in Gastroesophageal Adenocarcinoma: Guideline from the College of American Pathologists, American Society for Clinical Pathology, and American Society of Clinical Oncology. *Arch Pathol Lab Med*. 2016;140(12):1345-1363.

****REVISED** 12/26/2024**

CYG.48399 Validation/Verification - Predictive Marker Testing

Phase II



Predictive marker testing by in situ hybridization is validated/verified and records of validation/verification are retained.

*NOTE: For validation of **laboratory-developed and modified FDA-cleared/approved HER2 (ERBB2) breast predictive assays**, the validation must be performed on a minimum of 40 cases (20 positive and 20 negative samples).*

*For verification of **unmodified FDA-cleared/approved HER2 (ERBB2) breast predictive assays**, the laboratory must follow the instructions provided by the manufacturer. If the instructions do not list a minimum number of samples for assay verification, the verification must be performed on a minimum of 20 positive and 20 negative tissues.*

*For **other predictive marker assays**, the laboratory director must determine the appropriate number of positive and negative samples to be used to adequately validate/verify the test. In general, laboratories should consider using higher numbers of test cases when assessing laboratory-developed tests or modified FDA-cleared/approved tests than is necessary for unmodified FDA-cleared/approved tests for the same analyte. For genetic abnormalities where positive cases are rare, the laboratory director may determine that fewer validation cases are necessary. However, the rationale for using fewer cases must be recorded.*