

*data developed by these programs can provide a useful benchmark against which laboratory performance can be evaluated.*

**Evidence of Compliance:**

- ✓ Records of enrollment/participation in an educational peer-comparison program for interpretive flow cytometry **OR** records for participation in a laboratory-developed program circulating cases with other laboratories or within the laboratory's own practice with records of peer review



**REFERENCES**

- 1) Clinical and Laboratory Standards Institute (CLSI). *Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline—Second Edition*. CLSI document H43-A2. Clinical and Laboratory Standards Institute, Wayne, PA; 2007.

## QUALITY MANAGEMENT

### REAGENTS

#### Inspector Instructions:

	<ul style="list-style-type: none"> <li>• Sampling of new antibody validation records</li> <li>• Sampling of new lot/shipment antibody and detection system reagent confirmation records</li> </ul>
	<ul style="list-style-type: none"> <li>• What procedure does your laboratory follow to ensure manufacturer's recommendations are followed regarding the use of kit reagents/controls?</li> <li>• How do you confirm the acceptability of new reagent lots?</li> <li>• How do you evaluate the performance of newly prepared antibody cocktails?</li> </ul>

*Additional requirements are in the REAGENTS section of the All Common Checklist. Reporting requirements for use of analyte-specific reagents and other reagents used in laboratory-developed tests are included in the All Common Checklist (COM.40850).*

#### FLO.23250 Reagent Usage

**Phase II**

**The laboratory follows manufacturer's instructions for the proper use of reagents and controls or provides validation records if alternative procedures are used.**

**Evidence of Compliance:**

- ✓ Records of method validation if alternative procedures are used

**REFERENCES**

- 1) Caldwell CW. Analyte-specific reagents in the flow cytometry laboratory. *Arch Pathol Lab Med*. 1998;122:861-864

#### FLO.23275 Antibody Validation

**Phase II**

**The laboratory has records of initial validation of new antibodies prior to use in patient diagnosis.**

*NOTE: Antibodies used are validated on the cell sub-population of interest in the context of the antibody combination used in an assay.*