

**High quality materials with test system and matrix-appropriate target values are used for calibration and calibration verification whenever possible.**

*NOTE: Calibration and calibration verification must have defined analyte target values and appropriate matrix characteristics for the clinical specimens and specific assay method. Many instrument systems require calibration materials with system-specific target values to produce accurate results for clinical specimens.*

*Suitable materials for calibration verification include, but are not limited to:*

1. *Calibrators used to calibrate the analytical system*
2. *Materials provided by the manufacturer for the purpose of calibration verification*
3. *Previously tested unaltered patient/client specimens*
4. *Primary or secondary standards or reference materials with matrix characteristics and target values appropriate for the method*
5. *Third party general purpose reference materials that are suitable for verification*

*In general, routine control materials and proficiency testing materials are not suitable for calibration verification, except in situations where the material has been shown to be suitable (eg, specifically designated by the method manufacturer) or no other materials are available.*

**Evidence of Compliance:**

- ✓ Records of calibration and calibration verification

**REFERENCES**

- 1) ISO 17511:2020 In vitro diagnostic medical devices, Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples. International Organization for Standardization. 2020.
- 2) Department of Health and Human Services, Centers for Medicare and Medicaid Services. Clinical laboratory improvement amendments of 1988; final rule. *Fed Register*. 2003(Jan 24): [42CFR493.1255]

**CBG.12100 Recalibration/Calibration Verification Criteria**

**Phase II**



**Criteria for the frequency and acceptability of recalibration or calibration verification are defined and followed.**

*NOTE: Laboratories must either recalibrate or perform calibration verification at least every six months and if any of the following occur:*

1. *At changes of reagent lots unless the laboratory can demonstrate that the use of different lots does not affect the accuracy of patient/client results*
2. *If QC shows an unusual trend or shift or is outside acceptable limits and the system cannot be corrected to bring control values into the acceptable range*
3. *After major preventive maintenance or change of a critical instrument component*
4. *When recommended by the manufacturer*

**Evidence of Compliance:**

- ✓ Records of calibration verification documented at defined frequency

**REFERENCES**

- 1) Department of Health and Human Services, Centers for Medicare and Medicaid Services. Clinical laboratory improvement amendments of 1988; final rule. *Fed Register*. 2003(Jan 24):3707[42CFR493.1255(b)(3)]
- 2) Miller WG. Quality control. In: Henry's Clinical Diagnostic and Management by Laboratory Methods, 21st Edition, ed McPherson RA, Pincus MR. Saunders Elsevier, 2007: 99-111.

**CBG.12200 Recalibration**

**Phase II**

**The test system is recalibrated when calibration verification fails to meet the established criteria of the laboratory.**

**Evidence of Compliance:**

- ✓ Records of recalibration, if calibration or calibration verification has failed

**REFERENCES**

- 1) Department of Health and Human Services, Centers for Medicare and Medicaid Services. Clinical laboratory improvement amendments of 1988; final rule. *Fed Register*. 2003(Jan 24): [42CFR493.1255(a)(3)]