

DISCOVER

- Review a sampling of QC data over the previous two-year period. Select several occurrences in which QC is out of range and follow records to determine if the steps taken follow the laboratory procedure or corrective action

IMM.33930 QC - Waived Tests**Phase II**

The laboratory follows manufacturer's instructions for quality control, reviews results, and records acceptability prior to reporting patient results.

NOTE: Quality control must be performed according to manufacturer instructions. To detect problems and evaluate trends, testing personnel or supervisory staff must review quality control data on days when controls are run prior to reporting patient results. The laboratory director or designee must review QC data at least monthly or more frequently if specified in the laboratory QC policy.

*With respect to internal controls, acceptable control results must be recorded at a minimum, once per day of patient testing for each device.**

**Acceptable internal control results need not be recorded, if (and only if) an unacceptable instrument control automatically locks the instrument and prevents release of patient results.*

Evidence of Compliance:

- ✓ Records showing confirmation of acceptable QC results

IMM.33940 QC Corrective Action - Waived Tests**Phase II**

The laboratory performs and records corrective action when control results exceed defined acceptability limits.

CONTROLS - NONWAIVED TESTS

Inspector Instructions:

READ

- Sampling of quality control policies and procedures
- Sampling of QC records, including external and internal quality control processes

ASK

- How do you determine when quality control is unacceptable and when corrective actions are needed?
- How does your laboratory verify or establish acceptable quality control ranges?
- How does your laboratory perform QC for test procedures that report results as reactive, weakly reactive and nonreactive?
- What is your course of action when you perform test procedures that do not have commercially available calibration or control materials?

DISCOVER

- Review a sampling of QC data over the previous two-year period. Select several occurrences in which QC is out of range and follow records to determine if the steps taken follow the laboratory policy for corrective action