




Inspector Instructions:

| | |
|---|--|
|  | <ul style="list-style-type: none"> • Sampling of abnormal hemoglobin policies and procedures • Sampling of patient reports (confirmatory testing, comments) • Sampling of QC records |
|  | <ul style="list-style-type: none"> • Hemoglobin separation patterns (appropriate separations and controls) • Examine a sampling of medium (media) used to identify hemoglobin variants including alkaline/acid electrophoresis, isoelectric focusing, HPLC or other method |
|  | <ul style="list-style-type: none"> • What is your course of action when the primary screening method appears to show Hb S? • What is your course of action when the primary Hb screening method shows Hb variants migrating in non-A/non-S positions? |

CBG.20165 Hb S Primary Screen

Phase II

For samples that have screening results demonstrating the presence of a hemoglobin consistent with Hb S and suggesting a possible clinically significant condition, reporting of the screening results includes a recommendation that confirmatory testing be performed.

NOTE: For primary definitive diagnosis screening by electrophoresis or other separation methods, all samples with hemoglobins migrating in the "S" positions or peak must be tested for sickling hemoglobin(s). Known sickling and non-sickling controls both must be included with each run of patient specimens tested.

REFERENCES

- 1) Fairbanks VF. Hemoglobinopathies and thalassemias. Laboratory methods and case studies. New York, NY: BC Decker, 1980
- 2) Beuzard Y, *et al.* Isoelectric focusing of human hemoglobins, In Hanash, Brewer, eds. Advances in hemoglobin analysis. New York, NY: Alan R. Liss, 1981:177-195
- 3) Cossu G, *et al.* Neonatal screening of betathalassemias by thin layer isoelectric focusing. *Am J Hematol.* 1982;13:149
- 4) Bunn HF, Forget BG. Hemoglobin: molecular, genetic and clinical aspects. Philadelphia, PA: WB Saunders, 1986
- 5) Honig GR, Adams JG III. Human hemoglobin genetics. Vienna, Austria: SpringerVerlag, 1986
- 6) Jacobs S, *et al.* Newborn screening for hemoglobin abnormalities. A comparison of methods. *Am J Clin Pathol.* 1986;85:713-715
- 7) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part I. The introduction and thalassemia syndromes. *Lab Med.* 1987;18:368-372
- 8) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part II. The sickle cell disorders. *Lab Med.* 1987;18:441-443
- 9) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part III. Nonsickling disorders and cord blood screening. *Lab Med.* 1987;18:513-518
- 10) Armbruster DA. Neonatal hemoglobinopathy screening. *Lab Med.* 1990;21:815-822
- 11) Adams JG III, Steinberg MH. Analysis of hemoglobins, In Hoffman R, *et al.*, eds. Hematology: basic principles and practice. New York, NY: Churchill Livingstone, 1991:1815-1827
- 12) Mallory PA, *et al.* Comparison of isoelectric focusing and cellulose acetate electrophoresis for hemoglobin separation. *Clin Lab Sci.* 1994;7:348- 352
- 13) Awalt E, *et al.* Tandem mass spectrometry (MS) – A screening tool for hemoglobinopathies. *Clin Chem.* 2001;47(suppl):A165
- 14) Bradley CA, Kelly A. Comparison of high performance liquid chromatography with electrophoresis for measurement of hemoglobins A, A2, S, F, and C. *Clin Chem.* 2001;47(suppl):A172
- 15) Bradley CA, Kelly A. Calibration verification of hemoglobins A, A2, S, and F with an automated chromatography system. *Clin Chem.* 2001;47(suppl):A17315

CBG.20170 Daily QC - Hgb Separation

Phase II

Controls containing at least three known major hemoglobins, including Hb F and both a sickling and a nonsickling hemoglobin (eg, A, F and S) are applied with the patient specimen(s) and separations are satisfactory.

NOTE: There are written procedures for instruments with multiple electrophoretic chambers or capillaries to ensure that QC is performed on each individual chamber or capillary.

Evidence of Compliance:

- ✓ QC records reflecting the use of appropriate controls **AND**
- ✓ Data, tracings, or other appropriate testing results (eg, photographs, gels, cards with banding patterns) demonstrate appropriate controls and separation

REFERENCES

- 1) Fairbanks VF. Hemoglobinopathies and thalassemias. Laboratory methods and case studies. New York, NY: BC Decker, 1980
- 2) Beuzard Y, *et al.* Isoelectric focusing of human hemoglobins. In Hanash, Brewer, eds. Advances in hemoglobin analysis. New York, NY: Alan R. Liss, 1981:177-195
- 3) Cossu G, *et al.* Neonatal screening of betathalassemias by thin layer isoelectric focusing. *Am J Hematol.* 1982;13:149
- 4) Bunn HF, Forget BG. Hemoglobin: molecular, genetic and clinical aspects. Philadelphia, PA: WB Saunders, 1986
- 5) Honig GR, Adams JG III. Human hemoglobin genetics. Vienna, Austria: SpringerVerlag, 1986
- 6) Jacobs S, *et al.* Newborn screening for hemoglobin abnormalities. A comparison of methods. *Am J Clin Pathol.* 1986;85:713-715
- 7) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part I. The introduction and thalassemia syndromes. *Lab Med.* 1987;18:368-372
- 8) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part II. The sickle cell disorders. *Lab Med.* 1987;18:441-443
- 9) Fishleder AJ, Hoffman GC. A practical approach to the detection of hemoglobinopathies: part III. Nonsickling disorders and cord blood screening. *Lab Med.* 1987;18:513-518
- 10) Armbruster DA. Neonatal hemoglobinopathy screening. *Lab Med.* 1990;21:815-822
- 11) Adams JG III, Steinberg MH. Analysis of hemoglobins. In Hoffman R, *et al.*, eds. Hematology: basic principles and practice. New York, NY: ChurchillLivingstone, 1991:1815-1827
- 12) Mallory PA, *et al.* Comparison of isoelectric focusing and cellulose acetate electrophoresis for hemoglobin separation. *Clin Lab Sci.* 1994;7:348-352
- 13) Awalt E, *et al.* Tandem mass spectrometry (MS) – A screening tool for hemoglobinopathies. *Clin Chem.* 2001;47(suppl):A165
- 14) Bradley CA, Kelly A. Comparison of high performance liquid chromatography with electrophoresis for measurement of hemoglobins A, A2, S, F, and C. *Clin Chem.* 2001;47(suppl):A172
- 15) Bradley CA, Kelly A. Calibration verification of hemoglobins A, A2, S, and F with an automated chromatography system. *Clin Chem.* 2001;47(suppl):A17315

CBG.20175 Hemoglobin Variants

Phase II



All samples with hemoglobin variants not appearing to be Hb A, Hb F, or Hb S by the separation procedure in use, are reported with a recommendation to obtain confirmatory testing consistent with local screening program recommendations.

NOTE: The laboratory must have a defined process for the reporting of abnormal hemoglobin variants, developed in consultation with hematology advisors. The procedure must address the confirmatory testing recommendations, if confirmatory testing is not performed on-site.

Evidence of Compliance:

- ✓ Patient reports and records reflecting adherence to laboratory reporting procedures

REFERENCES

- 1) Giordino PC. Strategies for basic laboratory diagnostics of the hemoglobinopathies in multi-ethnic societies: interpretation of results and pitfalls. *Int J Hematol.* 2013;35:465-79.
- 2) Sabath DE. Molecular diagnosis of thalassemias and hemoglobinopathies: an ACLPS critical review. *Am J Clin Pathol.* 2017;148:6-15.
- 3) Greene DN, Vaughn CP, Crews BO, Agarwal AM. Advances in detection of hemoglobinopathies. *Clinica Chimica Acta.* 2015;15:439-50.
- 4) Troxler H, Kleinert P, Schmugge M, Speer O. Advances in hemoglobinopathy detection and identification. *Adv Clin Chem.* 2012;57:1-28.
- 5) Bain BJ. Haemoglobinopathy diagnosis: algorithms, lessons and pitfalls. *Blood Rev.* 2011;25(5):205-13.
- 6) Awalt E, *et al.* Tandem mass spectrometry (MS) – A screening tool for hemoglobinopathies. *Clin Chem.* 2001;47(suppl):A165
- 7) Bradley CA, Kelly A. Comparison of high performance liquid chromatography with electrophoresis for measurement of hemoglobins A, A2, S, F, and C. *Clin Chem.* 2001;47(suppl):A172
- 8) Bradley CA, Kelly A. Calibration verification of hemoglobins A, A2, S, and F with an automated chromatography system. *Clin Chem.* 2001;47(suppl):A173